

## DETERMINATION OF STRENGTHS OF AN ACID USING PH METER

Aim:-

To find out the strength of given hydrochloric acid solution by titrating it against sodium hydroxide (0.1N) using pH meter.

Principle -

When an alkali is added to an acid soln, the pH of the solution increases slowly, but at vicinity of the end point the rate of change of pH of the solution is very rapid. From the sharp break in the curve, we can find out the end point, from which the strengths of HCl can be calculated.

PROCEDURE -

- First standardize the pH meter using different buffer of known pH. Wash the glass electrode and reference electrode with distilled water and then with the acid soln.

- ② The given acid is diluted to 100mL using distilled water. 10mL of this made up solution is pipetted out into a 250mL beaker and 100mL of distilled water is added to it, so that the glass electrode and reference electrode are completely dipped.
- ③ Note the initial pH of the pure acid solution. Fill the burette with standard NaOH solution and run down into the beaker in small increments (0.5mL).
- ④ Stir the solution well using glass rod. Note down the pH after every successive addition. Continue the titration till beyond neutralisation point as indicated by an abrupt change in pH (at least in 3D int.).
- ⑤ Plot a pilot graph of volume of NaOH vs pH. The midpoints of the 's' shaped curve of the graph gives the equivalence point.
- ⑥ In order to get an accurate end point, perform one more similar titration.

Experiment No. \_\_\_\_\_ Name : \_\_\_\_\_

in a small volume range (1.0 mL on either side of the abrupt change in pH) and measure the pH after addition of every 0.2 mL of standard NaOH soln.

- ⑦. Plot a pair plot graph of volume of NaOH vs.  $\Delta \text{pH}/\text{mL}$ . Find out the exact end point from the pair plot. The break point of the pair plot graph gives us the end point.

### CALCULATION

Volume of NaOH,  $V_1 = 8.2 \text{ mL}$  (from FAIR graph)

Molarity of NaOH,  $N_1 = 0.1 \text{ N}$

Volume of HCl = 10 mL

Molarity of HCl,  $N_2 = 0.0485 \text{ N}$

$$\text{Then strength of HCl (N}_2\text{)} = \frac{N_1 \times V_1}{V_2} = \frac{0.1 \times 8.2}{10} = 0.082 \text{ N}$$

~~0.082 N~~


FAIR TITRATION (HCl vs NaOH).


S.	No.	Volume of NaOH ( $v$ ) (mL)	pH	SPN	$\Delta V$	$\Delta \text{SPN}/\Delta V$
1.		0	2.77	0	0	0
2.		7	3.60	0.83	0.2	4.15
3.		7.2	3.68	0.08	0.2	0.4
4.		7.4	3.83	0.15	0.2	0.75
5.		7.6	3.99	0.16	0.2	0.8
6.		7.8	4.27	0.28	0.2	1.4
7.		8.0	4.83	0.56	0.2	2.8
8.		8.2	5.80	0.97	0.2	4.85
9.		8.4	6.20	0.40	0.2	2.0
10.		8.6	6.56	0.36	0.2	1.8
11.		8.8	6.76	0.20	0.2	1
12.		9.0	6.95	0.19	0.2	0.95
13.		9.2	7.16	0.21	0.2	1.05
14.		9.4	7.44	0.28	0.2	1.4
15.		9.6	7.93	0.49	0.2	2.45
16.		9.8	8.69	0.76	0.2	3.8
17.		10.0	9.23	0.54	0.2	2.7
18.		10.2	9.82	0.59	0.2	2.95
19.		10.4	9.95	0.13	0.2	0.65
20.		10.6	10.05	0.10	0.2	0.5
21.		10.8	10.13	0.03	0.2	0.15
22.		11.0	10.24	0.11	0.2	0.55.

PILOT TITRATION. (HCl vs NaOH)

S <sub>2</sub> N <sub>o.</sub>	Volume of NaOH (v) (mL)	pH
1	0	2.85
2	1.0	2.88
3	2.0	2.90
4	3.0	2.94
5	4.0	3.01
6	5.0	3.10
7	6.0	3.22
8	7.0	3.41
9	8.0	3.80
10	9.0	6.03
11	10.0	6.75
12	11.0	7.66
13	12.0	9.69
14	13.0	10.11
15	14.0	10.42
16.	15.0	10.66

PILOT GRAPH.





Volume of NaOH (ml) →

RESULT →

The strength of given hydrochloric acid is 0.082 N.