

CHEMISTRY PRACTICAL

EXPERIMENT-I

DETERMINATION OF THE AMOUNT OF SODIUM CARBONATE AND SODIUM HYDROXIDE IN A MIXTURE BY TITRATION

Titration I: Standardisation of HCl

| S.No | Volume of Na_2CO_3 (mL) | Burette Reading (mL) | | Concordant Value (mL) | Indicator |
|------|---|----------------------|-------|-----------------------|---------------|
| | | Initial | Final | | |
| 1. | 20 | 0 | 18.7 | 18.7 | Methyl Orange |

Calculation:

$$\text{Volume of HCl} (V_2) = 18.7 \text{ mL}$$

$$\text{Normality of HCl} (N_2) = ? \text{ N}$$

$$\text{Volume of } \text{Na}_2\text{CO}_3 (V_1) = 20 \text{ mL}$$

$$\text{Normality of } \text{Na}_2\text{CO}_3 (N_1) = 0.05 \text{ N.}$$

$$N_1 V_1 = N_2 V_2$$



$$N_2 = \frac{N_1 V_1}{V_2}$$

$$= \frac{0.05 \times 20}{18.7}$$

$$N_2 = 0.05347$$

Titration II: Estimation of mixture ($\text{Na}_2\text{CO}_3 + \text{NaOH}$)

| S.No | Volume of the mixture (mL) | Burette Reading (mL) | | | Concordant Value (mL) | | Indicator |
|------|----------------------------|----------------------|-------------|------------|-----------------------|------|-----------------|
| | | Initial | Final - HPh | Final - Mo | HPh | Mo | |
| 1. | 20 | 0 | 23.5 | 32.5 | 23.5 | 32.5 | Phenolphthalein |
| 2. | 20 | 0 | 23.5 | 32.5 | | | Methyl Orange |

- A- Phenolphthalein end point- volume of acid used upto HPh end point.
 B- Methyl orange end point- Total volume of acid used till Mo end point.
 $(B-A) = C$ which is the Volume of acid used for $\frac{1}{2}\text{CO}_3^{2-}$ neutralization.
 $\therefore 2C =$ Total volume of acid used for complete neutralization of CO_3^{2-} ions.

Calculations:

1. Estimation of Na_2CO_3

$$\text{Volume of HCl (V}_1\text{)} = 2C \text{ mL} \\ = 18 \text{ mL}$$

$$\text{Normality of HCl (N}_1\text{)} = 0.05347 \text{ N}$$

$$C = B - A \\ = 32.5 - 23.5 \\ = 9 \\ 2C = 18 \text{ mL}$$

$$\text{Volume of the mixture (V}_2\text{)} = 20 \text{ mL}$$

$$\text{Normality of the mixture (N}_2\text{)} = ?$$

$$\text{N}_1\text{V}_1 = \text{N}_2\text{V}_2 \\ \text{HCl} \qquad \text{Mixture}$$

$$\text{N}_2 = \frac{\text{N}_1\text{V}_1}{\text{V}_2} = \frac{0.05347 \times 18}{20} \\ = 0.04812 \text{ N}$$

$$\therefore \text{Amount of } \text{Na}_2\text{CO}_3 \text{ present in the whole of the given solution} = \frac{\text{Normality of } \text{Na}_2\text{CO}_3 \times \text{Eq. wt of } \text{Na}_2\text{CO}_3}{10}$$

$$= \frac{0.04812 \times 53}{10}$$

$$= 0.2550 \text{ g/L}$$

2. Estimation of NaOH:

$$\begin{aligned}\text{Volume of HCl (V}_1\text{)} &= A - C \\ &= 23.5 - 9 \\ &= 14.5 \text{ mL}\end{aligned}$$

$$\text{Normality of HCl (N}_1\text{)} = 0.05347 \text{ N}$$

$$\text{Volume of mixture (V}_2\text{)} = 20 \text{ mL}$$

$$\text{Normality of mixture (N}_2\text{)} = ?$$

$$N_1 V_1 = N_2 V_2$$

HCl Mixture

$$N_2 = \frac{N_1 V_1}{V_2} = \frac{0.05347 \times 14.5}{20}$$
$$= 0.03876$$

$$\therefore \text{Amount of NaOH present in whole of the given solution} = \frac{\text{Normality of NaOH} \times \text{Eq wt of NaOH}}{10}$$
$$= \frac{0.03876 \times 40}{10}$$
$$= 0.1550 \text{ g/L}$$

RESULT:

Amount of Na_2CO_3 present in 100 mL of the given solution is 0.2550 g/L.

Amount of NaOH present in 100 mL of the given solution is 0.1550 g/L.

EXPERIMENT - II

ESTIMATION OF HARDNESS OF A WATER SAMPLE BY EDTA METHOD

Titration 1: Standardisation of EDTA

| S.No | Volume of Std. Hard Water (mL) | Burette Reading (mL) | | Volume of EDTA (V_1) (mL) | Indicator |
|------|-----------------------------------|----------------------|-------|----------------------------------|-----------|
| | | Initial | Final | | |
| 1. | 20 | 0 | 19.8 | 19.8 | EBT |

Calculation:

1 mL of Std. Hard water Contains 1 mg of CaCO_3 equivalent hardness.

∴ 20 mL of Std. Hard water Contains 20 mg of CaCO_3 equivalent hardness.

20 mL of Std hard water Consumes = V_1 mL of EDTA = 19.8 mL

∴

(V₁) 19.8 mL of EDTA = 20 mg of CaCO_3 equivalent Hardness.

1 mL of EDTA = $\frac{20}{19.8(V_1)}$ mg of CaCO_3 equivalent Hardness.

= 1.010 mg of CaCO_3 equivalent hardness.

Titration 2: Determination of Total Hardness

| S.No | Volume of Sample hard water (mL) | Burette Reading (mL) | | Volume of EDTA (V_2) (mL) | Indicator |
|------|--|----------------------|-------|-------------------------------------|-----------|
| | | Initial | Final | | |
| 1. | 20 | 0 | 12 | 12 | EBT |

Calculation:

20 mL of Sample hard water Consumes = V_2 mL of EDTA

= $V_2 \times 1$ mL of EDTA

= $V_2 \times \frac{20}{V_1}$ mg of CaCO_3 eqv. hardness

$\therefore 1000 \text{ mL of Sample hard water will consume} = V_2 \times \frac{20}{V_1} \times \frac{1000}{20} \text{ mg of CaCO}_3 \text{ eq. hardness}$

$$\begin{aligned}\therefore \text{Total hardness} &= 1000 \times \frac{V_2}{V_1} \text{ mg of CaCO}_3 \text{ equivalent hardness} \\ &= 1000 \times \frac{12}{19.8} \\ &= 606.0606 \text{ mg of CaCO}_3 \text{ equivalent hardness}\end{aligned}$$

Titration 3: Determination of Permanent hardness

| S.No | Volume of Boiled hard water (mL) | Burette Reading (mL) | | Volume of EDTA (V_3) (mL) | Indicator |
|------|----------------------------------|----------------------|-------|-------------------------------|-----------|
| | | Initial | Final | | |
| 1. | 20 | 0 | 8.5 | 8.5 | EBT |

Calculations:

$$\begin{aligned}20 \text{ mL of boiled water Sample Consumes} &= V_3 \text{ mL of EDTA} \\ &= V_3 \times 1 \text{ mL of EDTA} \\ &= V_3 \times \frac{20}{V_1} \text{ mg of CaCO}_3 \text{ eq. hardness}\end{aligned}$$

$$\therefore 1000 \text{ mL of boiled water Sample Consumes} = V_3 \times \frac{20}{V_1} \times \frac{1000}{20} \text{ mg of CaCO}_3 \text{ eq. hardness}$$

$$\begin{aligned}\therefore \text{Permanent hardness} &= \frac{V_3}{V_1} \times 1000 \text{ mg of CaCO}_3 \text{ eq. hardness} \\ &= \frac{8.5}{19.8} \times 100 \\ &= 429.2929 \text{ mg of CaCO}_3 \text{ equivalent hardness.}\end{aligned}$$

$$\text{Total hardness} = \text{Temporary hardness} + \text{Permanent hardness}$$

$$\begin{aligned}\therefore \text{Temporary hardness} &= \text{Total hardness} - \text{Permanent hardness} \\ &= 606.0606 - 429.2929 \\ &= 176.7677 \text{ mg of CaCO}_3 \text{ equivalent hardness.}\end{aligned}$$

RESULT:


Total hardness = 606.0606 ppm, Permanent hardness = 429.2929 ppm and Temporary hardness = 176.7677 ppm.

EXPERIMENT- IIIDETERMINATION OF STRENGTH OF AN ACID BY CONDUCTOMETRY

Titration 1: Titration between Std. NaOH and unknown HCl (Pilot)

| S.No | Volume of NaOH (mL) | Conductance (m ohm^{-1}) |
|------|---------------------|-------------------------------------|
| 1. | 0 | 4.5 |
| 2. | 1 | 3.8 |
| 3. | 2 | 3.2 |
| 4. | 3 | 2.7 |
| 5. | 4 | 2.3 |
| 6. | 5 | 2.0 |
| 7. | 6 | 1.91 |
| 8. | 7 | 1.83 |
| 9. | 8 | 1.69 |
| 10. | 9 | 1.52 |
| 11. | 10 | 1.41 |
| 12. | 11 | 1.56 |
| 13. | 12 | 1.68 |
| 14. | 13 | 1.79 |
| 15. | 14 | 1.84 |
| 16. | 15 | 1.92 |

GRAPH I: Titration between Std. NaOH and unknown HCl (Pilot)



Titration 2: Titration between Std. NaOH and unknown HCl (Fair)

| S.No | Volume of NaOH (mL) | Conductance (m ohm ⁻¹) |
|------|---------------------|------------------------------------|
| 1. | 8 | 1.72 |
| 2. | 8.2 | 1.70 |
| 3. | 8.4 | 1.69 |
| 4. | 8.6 | 1.67 |
| 5. | 8.8 | 1.63 |
| 6. | 9 | 1.50 |
| 7. | 9.2 | 1.48 |
| 8. | 9.4 | 1.45 |
| 9. | 9.6 | 1.42 |
| 10. | 9.8 | 1.39 |
| 11. | 10 | 1.35 |
| 12. | 10.2 | 1.31 |
| 13. | 10.4 | 1.28 |
| 14. | 10.6 | 1.36 |
| 15. | 10.8 | 1.45 |
| 16. | 11 | 1.49 |
| 17. | 11.2 | 1.54 |
| 18. | 11.4 | 1.58 |
| 19. | 11.6 | 1.63 |
| 20. | 11.8 | 1.67 |
| 21. | 12 | 1.70 |

Calculations:

$$\text{Volume of HCl (V₁) = 10 mL}$$

$$\text{Normality of HCl (N₁) = ?}$$

$$\text{Volume of NaOH (V₂) = 10.4}$$

$$\text{Normality of NaOH (N₂) = 0.1 N}$$

$$\frac{N_1 V_1}{\text{HCl}} = \frac{N_2 V_2}{\text{NaOH}}$$


$$N_1 = \frac{N_2 V_2}{V_1} = \frac{0.1 \times 10.4}{10}$$

$$N_1 = 0.104 N$$

RESULT:

The strength of the given HCl solution = 0.104 N

GRAPH 2: Titration between Std. NaOH and unknown HCl (Fair)




EXPERIMENT - IV

ESTIMATION OF THE AMOUNT OF Fe^{2+} ION USING $\text{K}_2\text{Cr}_2\text{O}_7$ BY POTENTIOMETRIC TITRATION

Titration 1: Titration between STD. $\text{K}_2\text{Cr}_2\text{O}_7$ and unknown FAS (Pilot)

| SNo | Volume of $\text{K}_2\text{Cr}_2\text{O}_7$ (ML) | EMF (mv) |
|-----|--|----------|
| 1. | 0 | 430 |
| 2. | 1 | 435 |
| 3. | 2 | 442 |
| 4. | 3 | 450 |
| 5. | 4 | 458 |
| 6. | 5 | 463 |
| 7. | 6 | 470 |
| 8. | 7 | 583 |
| 9. | 8 | 600 |
| 10. | 9 | 612 |
| 11. | 10 | 618 |
| 12. | 11 | 623 |
| 13. | 12 | 627 |
| 14. | 13 | 630 |
| 15. | 14 | 634 |
| 16. | 15 | 638 |


GRAPH 1: Titration between std $K_2Cr_2O_7$ and Unknown FAS (Pilot)



Titration 2: Titration between Std. $K_2Cr_2O_7$ and unknown FAS (Fair)

| S.N | Vol. of $K_2Cr_2O_7$ (mL) | EMF (mV) | ΔE | ΔV | $\Delta E / \Delta V$ |
|-----|------------------------------|----------|------------|------------|-----------------------|
| 1. | 4 | 450 | - | - | - |
| 2. | 4.2 | 452 | 2 | 0.2 | 10 |
| 3. | 4.4 | 456 | 4 | 0.2 | 20 |
| 4. | 4.6 | 460 | 4 | 0.2 | 20 |
| 5. | 4.8 | 463 | 3 | 0.2 | 15 |
| 6. | 5 | 468 | 5 | 0.2 | 25 |
| 7. | 5.2 | 470 | 2 | 0.2 | 10 |
| 8. | 5.4 | 473 | 3 | 0.2 | 15 |
| 9. | 5.6 | 478 | 5 | 0.2 | 25 |
| 10. | 5.8 | 481 | 3 | 0.2 | 15 |
| 11. | 6.0 | 485 | 4 | 0.2 | 10 |
| 12. | 6.2 | 489 | 4 | 0.2 | 10 |
| 13. | 6.4 | 493 | 4 | 0.2 | 10 |
| 14. | 6.6 | 497 | 4 | 0.2 | 10 |
| 15. | 6.8 | 500 | 3 | 0.2 | 15 |
| 16. | 7 | 573 | 73 | 0.2 | 365 |
| 17. | 7.2 | 584 | 11 | 0.2 | 10.8 |
| 18. | 7.4 | 589 | 5 | 0.2 | 25 |
| 19. | 7.6 | 591 | 2 | 0.2 | 10 |
| 20. | 7.8 | 593 | 2 | 0.2 | 10 |
| 21. | 8 | 596 | 3 | 0.2 | 15 |

GRAPH 2: Titration between STD. $K_2Cr_2O_7$ and unknown FAS (Fair)



Calculations:

Volume of FAS (V_1) = 10 mL

Normality of FAS (N_1) = ? N

Volume of $K_2 Cr_2 O_7$ (V_2) =

Normality of $K_2 Cr_2 O_7$ (N_2) = 0.1 N

$$N_1 V_1 = N_2 V_2$$

$$N_1 = \frac{N_2 V_2}{V_1}$$

$$= \frac{0.1 \times 7}{10}$$

$$N_1 = 0.07$$

$$\therefore \text{Amount of } Fe^{2+} \text{ present in the whole of the given solution} = \frac{\text{Normality of } Fe^{2+} \times \text{Eq. weight of } Fe^{2+}}{10}$$
$$= \frac{0.07 \times 55.85}{10}$$
$$= 0.03909 \text{ g/L}$$

RESULT:

The amount of Fe^{2+} present in the whole of the given solution = 0.03909 g/L